1. Product Name

Mirrored Fibre Optic Panel by ROMPA®

2. Product Code

18539

3. Colour

Black or white

4. Brief Description

100 strands of vinyl covered fibre optic light which gently change colour. Safe to handle and ideal for people with sight impairment as the optic fibre with smooth crimped ends can be held close to the eye.

5. Contents

- 1 x Panel

6. Snoezelen® Stimulations

- Sight

7. Best Use

Use in a darkened environment, for maximum visual impact. For maximum wireless interactivity, this panel can be used as part of ROMPA®'s Wi Fi system or as part of a room with ROMPA®'s Sensory Magic Touch Screen.

8. Compatible ROMPA® Products

19847 – 8 Colour Wirefree Controller

19939 – Wi Fi Interactive Carpet switch Set by ROMPA®

18870 - Talking Cube by ROMPA®

20402 – Mini Talking Cube

20563 – Wi Fi Microphone by ROMPA®

19214 – Colour Changing Panel by ROMPA®

18518 – Single Porter

18522 – Double Porter

9. Starting Up

- Remove the protective film
- Plug into mains socket. Keep mains cable out of reach

10. Detailed Description

100 strands of vinyl covered fibre optic light which gently change colour. Safe to handle and ideal for people with sight impairment as the optic fibre with smooth crimped ends can be held close to the eye. Set in a panel against an acrylic mirrored backing. The frame is made of textured plastic.

11. Safety


- Always supervise the use of this product.
- For indoor use only.
- Fibre optics should never be placed in the mouth.

12. Technical Specification

Size: 115 x 72.5 x 11 cm
Weight: 21kg
Input: 100-240v AC 50/60Hz 0.25A
Output: 12v DC 1.25A 15W max
Cover: phthalate-free
CE: This panel meets the requirements of BS EN 9001-2000

13. Installation

1. Choose a safe and suitable location for your panel. The panel is ONLY to be used in dry environments.
2. No fixings are included in this product for the wall. When fixing it to the wall please consult a suitably qualified professional builder or from an estates department to understand the construction of the wall the panel is to be fixed to, so that the correct fittings can be used.
3. Attach fixings such as mirror plates to the back of the panel. Check the fixings e.g. mirror plates you buy can easily support 21kg. You may wish to attach two mirror plates to the back of the panel at the top and two on the back at the bottom. When attaching fixings to the back of the panel please proceed with caution – be aware that there are electrical assemblies inside the panel.

4. Alternatively, this Panel may be used in a Porter 18518/18522 – see ROMPA® catalogue/website for further information.
5. Remove the clear protective film from the panel front.

6. Plug in at the wall.

14. Care and Maintenance

- Wipe the exterior of the Panel with a clean, soft and very slightly dampened cloth. Wipe the fibre optic strands with a soft, damp cloth. Do not allow the Panel to become wet. Do not immerse the fibre optics in water.
- Do not strike the Panel – handle it with appropriate care.
- There should be no need to access the interior of the panel. The light source inside the panel contains LEDs which last thousands of hours.
- In the highly unlikely event of product failure, please contact ROMPA®.

15. Troubleshooting

My Panel does not light up?

- Ensure that the panel is plugged in
- Check that the mains power lead connector is pushed firmly into the light source.
- Check the fuse that is in the Mains Power connector fitted to the back of the light source. The fuse is in a small drawer that can be prised out only when the Mains Lead is removed. If the fuse needs replacing, it should be replaced with a 1Amp Time Delay Fuse. Dimensions 5mm x 20mm
- Check the fuse in the plug (if UK plug fitted). This should be a 5A standard mains fuse.
- Turn the light source OFF, wait 1 minute and then turn ON again

My Panel sometimes changes colour randomly, and then goes back to its pre-programmed sequence

- The panel may be picking up Wi Fi signals from somebody nearby who has a Wi Fi controller e.g. Colour Match Panel, 8 Colour Wirefree Controller etc.
- If you have 2 neighbouring rooms each using Wi Fi products, then the devices in both rooms can be controlled by any nearby Wi Fi controllers. An "All-In-One" Controller would be required to alter the wireless settings to operate without interfering with one another.

Insufficient Brightness from the Light Source

- Ensure that the Fibre Optic harness is pushed fully into the front of the light source. It will stop when it is fully pushed in. Tighten the top screw on the light port on the front of the light source to prevent the harness from being pulled out.

Further copies of these instructions can be downloaded at www.rompa.com

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at producthelp@rompa.com

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