1. Product Name

Swing Frame

2. Product Code

16884

3. Colour

Green

4. Brief Description

This swing frame is for indoor use and can accommodate toddlers, adults and youths with different sized swings (to be purchased separately)

5. Contents

1.	Base Rail	8	Hook
2	Leg	9	Button Head Bolt & Screw
3	Top Rail	10	End Fitting Bolt
4	End fitting	11	Screw & Nut for Hook
5	T-Clamp	12	Screw & Nut for Bracket
6	Bracket – upper piece	13	End Cap
7	Bracket – lower piece	14	2 x allen keys

Parts 1, 2, 3 and
4

6. Snoezelen® Stimulations

- Vestibular

7. Best Use

Use indoors to provide vestibular input for clients.

8. Compatible ROMPA® Products

Adult swing seat: 16881
Youth swing seat: 16882
Toddler swing seat: 16883

9. Starting Up

1. Unpack the product and dispose carefully of packaging materials.
2. Familiarise yourself with all the components shown and listed in the contents section.
3. Decide where the Swing Frame will stand – you will not want to move it far once it is assembled.
4. Position in a safe and appropriate place. Ensure that walkways and doors are not obstructed and that the item is not positioned near to sharp edges such as furniture corners etc to minimize potential risk if client/s move around.
5. Assemble item safely (see installation).
6. Decide if the client needs assistance transferring to use the item.

10. Detailed Description

This sturdy and reliable swing frame can accommodate one swing at a time. Three different sized swing seats can be attached making it inclusive for a range of people (see compatible ROMPA® products).

11. Safety

- **Always supervise use of this product.**
- This item is not a toy, it is a piece of sensory equipment. Do not allow unsupervised use.
- For indoor use only. Do not leave the item outdoors.
- For 3 years and over.
- The Swing Frame should not be mouthed.
- Remove all packaging before use and store or dispose of this appropriately.
- Keep and use away from sources of heat e.g. fire, radiators etc.
- Before assembly, consider the surface under and around the swing – you may want to consider matting, softer surfaces etc. Consider the use area and landing surface to the front, rear and sides of the swing – a minimum distance of twice the height of the Top Rail above the floor on all sides.
- There should be at least six feet of obstacle-free, safe area to the sides of the Swing Frame.
- The Swing Frame should be installed on a flat, stable surface, away from sharp, pointed or hazardous objects.
- Take time and care to assemble and install this product correctly. Failure to do so could result in injury.


This product is not a toy. Any activity involving motion, height, or rotation may cause serious injury. This product is intended for use under the supervision of trained professionals.

12. Technical Specification


Weight capacity: 500 lbs
Size: Approximately 230cmH x 250cmL x 250cmW
Material: Metal

13. Installation


- Attach top rail and legs to end fitting (part 4).
- Attach frame to base rails with T-clamps (part 5).
- Stand upright. Plumb and level the whole swing.
- Tighten all nuts and bolts. Attach swing hangers (parts 6 and 7).
- Locate Top Rail bolting. The Top Rail (part 3) must be bolted into each of the End Fittings (part 4).

A

B

14. Care and Maintenance

- Set up a schedule for checking the Swing Frame on a regular and frequent basis.
- Make sure it is adequately maintained, repairing or replacing any worn parts as appropriate.
- Establish a regular maintenance schedule to check the Swing Frame for signs of wear.
- Check MONTHLY for loose bolts, damaged or broken parts.
- Tighten and/or replace immediately.
- Replace any parts that show a reduction of 25% of original size, if not before this amount of wear.

15. Troubleshooting

Not applicable

Further copies of these instructions can be downloaded at www.rompa.com

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at producthelp@rompa.com

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