## VIBRATION PILLOW 16700-1, 16776-8, 16916, 16981

**NOTE:-** The Vibration Pillow may help ease aches and pains. However, it is not a cure or alternative to professional medical attention. If pain persists consult a doctor.  
Uses 2 x D batteries (not supplied)

### TO FIT BATTERIES

1. Open the zip on the back of the pillow.
2. Remove the plastic holder and insert the batteries observing the polarity indicated. Ensure that they are placed in the same direction.
3. Replace the holder and close the zip.

### USE OF PILLOW

Gentle pressure activates the sensor and causes the pillow to vibrate. The pillow can help soothe away aches and pains and increase blood circulation.

### SUGGESTED USES

The Vibration Pillow may help the following:

1. For gentle relaxation and to ease swollen ankles place feet on the pillow.
2. For soothing tired and aching muscles, place the pillow behind the calves.
3. For easing the discomfort of aches and pains in the joints, place the pillow behind the knees.
4. To ease sciatic pain, place on the lower back.
5. To relieve stiff necks, as well as tension and stress which can cause headaches, place the pillow between the shoulder blades.
6. To activate vibrations throughout the back region, place the pillow in the middle of the back.

### CARE & MAINTENANCE

Replace the batteries as required. The cover may be removed with care and hand-washed in cool water using a very mild cleaning product. Do not immerse in water.

### SAFETY

Always supervise the use of this product.

## SWITCHABLE VIBRATION PILLOW 15800, 18944

### OPERATING INSTRUCTIONS

1. Insert any ROMPA® non-electronic switch (not included) into the jack socket on the pillow.  
**NOTE:-** The pillow will no longer vibrate when pressure is applied.
2. Operate switch to cause the pillow to vibrate.
3. Remove the jack plug of the switch from the socket and normal operation can be resumed.

### CARE & MAINTENANCE

Replace the batteries as required. The cover of the Switchable Pillow cannot be removed – wipe clean only. Do not immerse in water.

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