GLITTER GEL SHAPES - SET OF 5 20092

CONTENTS

5 x Glitter Gel Shapes

BEFORE USE

1. You may wish to put a protective cloth on your table just in case of accidental damage to the shapes. Likewise, if not using a table, you may to cover the client's lap with a protective cloth just in case.
2. Remove all buckles, keys and any sharp objects before using so as to avoid accidental damage.

USE

Use under supervision.

Handle with appropriate care.

Suggestions for use

1. Use in colour and shape recognition games.
2. Use in pre-handwriting activities – encourage the client to isolate the index finger and trace the shape through the gel.
3. Gently move the gel around using the fingers.
4. Provide further tactile stimulation by changing the temperature of the shapes. To heat them, soak them in warm-hot water (not boiling); to cool the shapes, place them in the fridge. Do not freeze the shapes.

SAFETY

1. Always supervise the use of this product.
2. Do not place in the mouth. Not suitable for those who mouth.
3. If heating/cooling the shapes, make sure the shapes are a comfortable temperature before handing them to clients.

CARE AND MAINTENANCE

1. Handle with care.
2. Do not jump on the shapes or tread on them.
3. Store in a safe place, away from extremes of temperature and direct sunlight.
4. Wipe clean with clean water and a soft cloth.

SPECIFICATION

Filled with non toxic gel

Latex-free, tear-resistant and antibacterial vinyl exterior

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