



### **1. Product Name**

TheraGym® Over the Moon Swing Set A

### **2. Product Code**

21581

### **3. Colour**

As shown


### **4. Brief Description**


Provide many opportunities to experience multi-sensory input. Each set includes a half moon-shaped support made of strong Baltic birch wood, a heavy-duty bar strap and several options to secure these to your suspension system.

## 5. Contents

- 1 x Vertical Bouncer Kit:
  - 2 x Single Vertical Bouncer units
  - 2 x Double Vertical Bouncer units
  - 2 x locking Safety Snaps
- 1 x Bouncing Swing:
  - 1 x half-moon wood bar
  - 1 x heavy duty green bar strap
  - 1 x swing seat with 2 industrial-strength straps
- 1 x Sling Swing:
  - 1 x Sling Swing
  - 2 x locking Safety Snaps
- 1 x Safety Snap  
*This set includes 5 Safety Snaps in total*







## 6. Snoezelen® Stimulations

- Vestibular

## 7. Best Use

Use from a safe and appropriate suspension system.

## 8. Compatible ROMPA® Products

21928 – Hanging Bracket

21559 – Safety Rotational Device

23011 (21570, 21590) – Safety Snap

See [www.rompa.com](http://www.rompa.com) for mats and a large range of Sensory Integration products

## 9. Starting Up

- There are no short-cuts when setting up and using Sensory Integration equipment. Please keep these instructions in a safe place and refer to them before each use.

## 10. Detailed Description

### **Over the Moon Set A**

Provide many opportunities to experience multi-sensory input. Each set includes a half moon-shaped support made of strong Baltic birch wood, a heavy-duty bar strap and several options to secure these to your suspension system. Working load: 136kg (300 lbs).

Set A includes a Vertical Bouncer Kit with 2 Single Bouncer units; 2 Double Vertical Bouncer units; and 2 Safety Snaps (5 Safety Snaps in total in the set).

### **Sling Swing**

Provides stimulation in a prone or seated position with its versatile aerial hammock suspension and quality construction. Safe for children and adults. Includes one locking Safety Snap. This product can be suspended from one suspension point or from two, such as from the Half-moon wood bar of the Over the Moon Set A. Weight capacity: 136kg (300 lbs).

Size: 124cm L x 76cmW (49”L x 30”W)



### **Bouncing Swing**

A popular swing for home and clinical use for proprioceptive stimulation and postural work. Using the Vertical Bouncer Kit attachments, it is ideal for providing multiple sensory opportunities. The swing hangs 122cm (4') below the point of attachment. Working load: 68kg (150 lbs).



### Sling Swing

1. Slide the two Suspension Straps (located at the top of the swing) onto the Safety Snap.
2. Connect the Safety Snap to your point of suspension and lock the Safety Snap closed.



### Bouncing Swing

1. Pinch together one section of the green bar strap and thread it through the centre hole at the top of the half-moon bar.




2. Pass one portion of the green bar strap loop through the loop on the other side of the half-moon bar, and pull tight.



3. Connect a locking Safety Snap to the open green bar strap loop, then hook the Safety Snap to your point of suspension.



4. Hook the Safety Snaps from the swing or accessory you wish to attach to the holes on each opposite end of the half-moon base.




### **Vertical Bouncer Kit**

1. Open one Safety Snap on the Vertical Bouncer Kit and unhook it to flatten out the bouncer.





2. Select the desired single or double bouncer attachments and disconnect the unwanted accessories from the Safety Snap.



3. Attach one of the Vertical Bouncer Kit Safety Snaps to the reinforced loop on the top of the swing.



4. Connect the Safety Snap on the other end of the Vertical Bouncer Kit to your secured end point of suspension.



## **14. Care and Maintenance**

1. All the parts of your suspended products should be checked on a regular basis (at least monthly, and more if the swing is used often). Use a maintenance sheet (such as the example one provided) to record these checks.
2. Check all the parts for signs of wear and tear before each use and replace these as appropriate. For example, check the Sling Swing for loose stitching.
3. Use and store this product indoors, in a safe place.
4. Dust or wipe the components clean with a damp cloth. Do not immerse in water. Allow to dry before use.

## **15. Troubleshooting**

*Not applicable*

Further copies of these instructions can be downloaded at [www.rompa.com](http://www.rompa.com)

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at [producthelp@rompa.com](mailto:producthelp@rompa.com)

