



1. Product Name

TheraGym® Over the Moon Swing Set B

2. Product Code

21583

3. Colour

As shown

4. Brief Description

Provide many opportunities to experience multi-sensory input. Each set includes a half moon-shaped support made of strong Baltic birch wood, a heavy-duty bar strap and several options to secure these to your suspension system.

5. Contents

- 1 x Half-moon wooden bar



- 1 x Heavy-duty green bar strap



- 1 x Trapeze Bar attachment



- 2 x reinforced nylon straps



- 2 x Acrobat Rings



- 2 x Safety Snaps



Please note: the reinforced nylon straps may be supplied already attached to the Trapeze Bar attachment

6. Snoezelen® Stimulations

- Vestibular

7. Best Use

Use from a safe and appropriate suspension system.

8. Compatible ROMPA® Products

21928 – Hanging Bracket

21559 – Safety Rotational Device

23011 (21570, 21590) – Safety Snap

See www.rompa.com for mats and a large range of Sensory Integration products

9. Starting Up

- There are no short-cuts when setting up and using Sensory Integration equipment. Please keep these instructions in a safe place and refer to them before each use.

10. Detailed Description

Over the Moon Set B

Provide many opportunities to experience multi-sensory input. Each set includes a half moon-shaped support made of strong Baltic birch wood, a heavy-duty bar strap and several options to secure these to your suspension system. Working load: 136kg (300 lbs).

Set B includes an 81cm Trapeze Bar and a set of 2 Acrobat Rings (either can be attached to the half-moon-shaped support with 91cm reinforced nylon straps) and 2 Safety Snaps to attach the straps to the wooden support.

11. Safety

- Always supervise the use of this product
- Always use the Rotational Device to help prevent wear of your swing and the suspension points.
- Do not exceed the working load: 136kg (300 lbs).
- Sensory Integration products are intended to be used by professionals.
- Use safety mats on the floor, as well as on the walls as necessary.
- Use a safety helmet, and any other appropriate safety measures if required for your client's needs.
- Use only steel safety snaps as aluminium safety snaps wear more quickly.
- For indoor use only.
- Check all hardware and stitching before each use to check these are intact.

This product is not a toy. Any activity involving motion, height, or rotation may cause serious injury. This product is intended for use under the supervision of trained professionals.

12. Technical Specification

Working load: 136kg (300 lbs)
Size: Trapeze Bar 81cmL
Reinforced nylon straps 91cmL

13. Installation

Allow plenty of unobstructed space around your swing. Suspended equipment needs a minimum ceiling height of 2.5m (8 feet). Suspended equipment should be positioned a clear distance of 1.8m (6 feet) from walls and obstacles.

1. Pinch together one section of the green bar strap and thread it through the centre hole at the top of the half-moon bar.




2. Pass one portion of the green bar strap loop through the loop on the other side of the half-moon bar and pull tight.



3. Connect a locking Safety Snap to the green bar strap loop, then hook the Safety Snap to your suspension system.



- Hook the Safety Snaps from the swing or accessory you wish to attach to the holes on each opposite end of the half-moon bar.



- To attach the Acrobat Rings to the straps, thread one end of the strap through the Acrobat Ring.



- Pull one end of the strap through the loop on the other end of the strap.



7. Pull tight so that the threaded loop secures the ring on one end.



8. Attach the open loop end of the strap with a Safety Snap to the half-moon bar at one of the end holes.



9. To attach the Trapeze Bar, open the loop on one end of the strap and slide it on to one end of the Trapeze Bar, sliding it towards the middle of the bar just past the hole at the bar's end.



10. With the threaded loop to the inside of the hole in the bar, pinch the open loop end of the strap and thread it through the hole in the bar.





11. Pull the strap all the way through the hole in the bar and then slide the threaded loop to the outside of the bar so that it touches, but can't go past, the hanging strap. Pull the hanging strap straight and attach to the half-moon bar with a Safety Snap. Repeat for the other strap on the other side of the bar.



Safety Snaps

Our Safety Snaps and the included Carabiner can be locked. Turn the locking nut to open the snap (to slide this nut down to allow the snap to be opened). Turn the nut in the opposite direction to slide it up and over the opening to lock it.





14. Care and Maintenance

1. All the parts of your suspended products should be checked on a regular basis (at least monthly, and more if the swing is used often). Use a maintenance sheet (such as the example one provided) to record these checks.
2. Check all the parts for signs of wear and tear before each use and replace these as appropriate. For example, check the ropes and straps for loose stitching.
3. Use and store this product indoors, in a safe place.
4. Dust or wipe the components clean with a damp cloth. Do not immerse in water. Allow to dry before use.

15. Troubleshooting

Not applicable

Further copies of these instructions can be downloaded at www.rompa.com

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at producthelp@rompa.com

