



### **1. Product Name**

Chatsworth Range by ROMPA® with vibration

### **2. Product Codes**

22520, 22554 – Chair with vibration

22522, 22556 – Chair with high back and vibration

22524, 22558 – Settee 2 Seat with vibration

22526, 22560 – Settee 2 Seat with high back and vibration

22528, 22562 – Left Hand 2 Seat Settee with vibration

22530, 22564 – Left Hand 2 Seat Settee with high back and vibration

22532, 22566 – Right Hand 2 Seat Settee with vibration

22534, 22568 – Right Hand 2 Seat Settee with high back and vibration

22536, 22570 – Corner Unit with vibration

22538, 22572 – Corner Unit with high back and vibration

22544, 22578 – Left Hand 3 Seat Settee with vibration

22546, 22580 – Left Hand 3 Seat Settee with high back and vibration

22540, 22574 – Right Hand 3 Seat Settee with vibration

22542, 22576 – Right Hand 3 Seat Settee with high back and vibration

22548, 22582 – Settee 3 Seat with vibration

22550, 22584 – Settee 3 Seat with high back and vibration

### **3. Colour**

Available in a range of colours, please contact us for the full range.

#### **4. Brief Description**

A range of soft furniture available in cheerful, bright colours.

#### **5. Contents**

- 1 x piece of Chatsworth furniture with vibration

#### **6. Snoezelen® Stimulations**

- Touch
- Sight

#### **7. Best Use**

Use multiple pieces together as a suite.

#### **8. Compatible ROMPA® Products**

See [www.rompa.com](http://www.rompa.com) for the full range. To make the furniture even more comfortable and inviting-looking, why not add a Foot Stool and some Chatsworth Scatter Cushions? Choose matching colours or contrasting shades for visual stimulation.

#### **9. Starting Up**

- Unpack and dispose of packaging safely and appropriately
- Connect the DC plug from the furniture to the DC lead from the power supply
- Connect IEC plug to the power supply
- Plug in
- Sit on the furniture to activate the vibration

#### **10. Detailed Description**

The ever-popular Chatsworth range is now better than ever:

- Ergonomically designed
- With a higher cushion for additional back support
- Solid foam without a rigid wooden frame
- Multi-sensory, with eye-catching and tactile fabric options
- Contemporary and aesthetically pleasing for a multitude of environments, including homes and waiting areas
- Softer, slimmer and more curvaceous styling
- Supportive without compromising on comfort

Available in a huge choice of colours. Choose from our standard PVC or 2-tone options - Faux Leather combined with Tactile Fabric. As these products have seams, they are not completely waterproof. Where there is the occasional risk of incontinence amongst furniture users, we recommend that PVC is chosen.

The Chatsworth range includes individual pieces of furniture as well as settees and corner units that can be placed together to maximise the space available. These modular pieces can be used in a variety of combinations and are perfect for waiting areas. For example, team a 3 Seat settee that has a single arm with a Corner Unit to fill a whole wall. A second settee with the opposite single arm could be added to the other side of the Corner Unit to fill 2 adjacent right-angled walls. With so many colours, sizes and sensory options to choose from, it's up to you!

Remember, L/H refers to the left hand side from the perspective of the seated user, not from the view of someone looking at the settee in front of them.

## 11. Safety



Please read all the following safety information before each use

### **ALWAYS SUPERVISE THE USE OF THIS PRODUCT**

- Use on a flat, stable surface, capable of holding the product's weight and that of the people sat on it.
- Avoid sharp or pointed objects which could puncture the product. Take care with zips on clothing etc. Remove shoes.
- Remove all packaging and either store this in a safe place out of reach of clients or dispose of this responsibly. Plastic bags and wrapping can be hazardous to children and vulnerable people.
- Chatsworth Chairs contain a foam filling. Do not attempt to access the interior part of this product. Never allow children or vulnerable people to access the interior part of this product.
- Should any seams or the zip become damaged, contact ROMPA® immediately to arrange repair or replacement.
- Keep away from sources of heat and fire.

## **12. Technical Specification**

Size: 75cmH (standard height)  
Material: PVC / Textured Fabric and Faux Leather Fabric  
Power Supply: input 100-240V~ 47-63Hz 1.4A max  
Output 12V--- 5.0A

## **13. Installation**

1. Unpack.
2. Place the furniture indoors, on a flat stable surface that is capable of holding the product's weight and that of the people sat on it.
3. Connect the DC plug from the furniture to the DC lead from the power supply.
4. Connect the IEC plug to the power supply.
5. Plug in at the wall.
6. Switch on at the wall (there is no ON/OFF switch on the product).
7. Sit on the furniture to activate the vibration. Applying pressure to the furniture activates the vibration. The vibration modules are located in the back rest and seat area (not in the cushion).
8. The intensity of the vibration is not adjustable.

## **14. Care and Maintenance**

- Organise a regular cleaning schedule.
- Treat this furniture with care e.g. do not eat whilst on the furniture to avoid spillages.
- Look after your furniture to maintain its appearance and condition.

### **FABRIC MAINTENANCE**

If you look after your soft cushioning it should give you many years of service.

But please observe the following:-

1. Fabric may wear, tear or "rub off" if subjected to sharp objects or continual rubbing by hard objects. Avoid shoes, buckles or sharp/hard objects.
2. From time to time lift/move furniture to clean all surfaces and the floor underneath.

### **Faux Leather**

#### **General Cleaning:**

Visible and loose dirt can be picked up with a vacuum cleaner – do not use a brush type vacuum. To get rid of any dust or marks on the surface, the Faux Leather MUST be wiped DAILY with a moist soft cotton cloth. Use a soapy microfiber or cotton cloth and gently rub stains in circular movements. For more stubborn marks a dilute bleach solution can be used – 5% bleach to 95% water to give 10,000ppm of available chlorine. Chlor-Clean, Haz Tabs and Tristel Fuse for Surfaces

are suitable for use on this fabric. Continued use and stronger solutions may damage the surface and cause the vinyl to crack.

#### Stain Removal:

Blot excess liquid as quickly as possible with a clean microfiber or 100% cotton white absorbent cloth or sponge. If necessary, use clean warm water. Clean the entire area where the spot occurred on the cushion, arm etc. Water may dissolve dirt particles and making cleaning easier. Spills that are still moist are the easiest to remove. Leave to air dry naturally. Do not use a hair dryer etc. to speed up the drying.

#### *Coffee, tea etc. spills:*

Use an absorbent white cloth or piece of paper. DO NOT RUB. Allow the stain to be absorbed and removed from the material surface. Then rub gently with white microfiber cloth or tissue paper in circular movements to remove any remaining dampness and residue.

#### *Ink, ballpoint pen, Denim:*

These are almost impossible to remove once set. Rub the stain by means of gentle circular movements with a microfiber cloth or tissue dampened with alcohol and water (1:1). After this, rub the surface with gentle circular movement with tissue paper to dry it. Continued use of alcohol-based solutions may damage the surface and cause cracking.

#### *Butter, oil, grease stains:*

Wipe excess butter, oil or grease off the surface with a clean microfibre cloth or 100% cotton white dry cloth. Use a soapy cotton cloth and gently rub the stain with circular movements. Use a clean damp cloth to wipe the area and dry the surface.

#### *Chewing gum:*

Place ice cubes in a plastic bag and apply on to the chewing gum. Wait several minutes then lift the ice bag and remove the chewing gum.

DO NOT USE CHEMICALS, SADDLE SOAP, CLEANING SOLVENTS, FURNITURE POLISH, OILS, VARNISH, ABRASIVE CLEANERS OR AMMONIA IN WATER AS THESE WILL DAMAGE THE FABRIC.

It is crucial that any staining is treated immediately. Delays in cleaning will give stains time to react with the plasticizer in the faux leather, making it more and more difficult to remove.

**Please note:** These are recommended or suggested methods of cleaning. The manufacturer is not responsible for damage incurred while cleaning. Always try the cleaning method in a hidden area first to check the results. Remember, if you are unsure seek further advice.

### **Textured Fabric**

*General Cleaning:* Regularly vacuum and wipe with a damp cloth to remove dust particles or marks on the surface. Most stains can be removed using the wipe clean method.

*Spillages:* Remove quickly with an absorbent dry cloth.

*Minor soiling:* Wipe with a clean damp microfiber cloth.

*Water-based stains:* Use mild liquid detergent with warm water. Avoid excessive amounts of water. Rinse thoroughly then dry with an absorbent microfibre cloth. Use cold water for stains produced by bodily fluids. All stains caused by bodily fluids should be removed immediately.

*Other stains:* Use proprietary products, following the instructions on the container carefully.

All soap residue must be rinsed out of the fabric as it will attract further stains.

Hypochlorite or Alcohol based solutions can be used in extreme circumstances. Use half a cup of household bleach with 5 litres of water – 10,000ppm.

Steam cleaning can be used in situ.

If necessary the face side of the fabric only could be ironed at a medium temperature.



DO NOT DRY CLEAN.



It is crucial that any staining is treated immediately.

**Please note:** These are recommended or suggested methods of cleaning. The manufacturer is not responsible for damage incurred while cleaning. Always try the cleaning method in a hidden area first to check the results. Remember, if you are unsure seek further advice.

## 15. Troubleshooting

In the unlikely event of product failure, please check the following:

**Problem: No vibration**

- Check that the furniture is plugged in at the wall
- Check that the product is switched on at the wall

Further copies of these instructions can be downloaded at [www.rompa.com](http://www.rompa.com)

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at [producthelp@rompa.com](mailto:producthelp@rompa.com)