SOUTHPAW FROG SWING 29037

Take time to familiarise yourself with the use and maintenance of this equipment BEFORE you start to use it. Keep these instructions in a safe place, and download more copies as required from www.rompa.com to ensure that everyone using this equipment is aware of these instructions.

IMPORTANT: please check equipment immediately for any damage/loss that may have occurred during transit.

NOTE: All hanging equipment requires a minimum ceiling height of 8 feet/2.5 metres and should be positioned a clear distance of 6 feet/1.8 metres from any adjacent walls and other obstacles.

Working load: 35-100 lb (150 lbs/68 kg with additional loops)

Suspension system not included - see ROMPA® Catalogue for details of Hanging Bracket (13672), Safety Snaps (13799), Safety Rotational Device (13738), Height Adjuster (13626) and Therapy Rope with Eye Splice (13627, 14535). It is your responsibility to ensure that your ceiling can support this equipment and use of it. The use of a Southpaw Height Adjuster 13626 is advisable. Use of a Southpaw Safety Rotational Device 13738 is essential to prevent damage to equipment and help ensure its safe use. Keep all parts together.

CONTENTS

1	x	Frog Swing
3	x	Stimulation Loops
2	x	Locking Safety Snaps

SUSPENSION

Attach your Frog Swing to your ceiling support system using the Locking safety Snap provided. The Frog Swing is now ready for use.

MAINTENANCE CHECK LIST

All of the following wear points should be checked on a monthly basis and recorded on the checklist attached.

- Locking Safety Snaps: see Safety Snap instructions
- Stimulation Loops: check for loose stitching around the 'ball' near the top and bottom, and holding the loops on the top and bottom. If found, discontinue use till repaired or replaced.
- White Poly Ropes: check for fraying rope. If found, discontinue use until repaired.
- Rope splice-clips: these are the small metal rings found near each splice of the white poly rope. If any are found to be loose or missing, discontinue use.
- Wood beam: check for cracks or splinters. If cracked, discontinue use immediately. For splinters, sand smooth and refinish with clear lacquer.

STORAGE

Store in a safe place out of reach so as to avoid unwanted access by clients.

Store out of direct sunlight.

Keep away from chemicals and household cleaning products.

Keep away from sources of extreme heat

IMPORTANT: KEEP ALL PARTS TOGETHER

CLEANING

Clean with a soft, slightly damp clean cloth

SAFETY

Always supervise the use of this product

For indoor use only, in dry environments.

Take great care to ensure that Southpaw products are appropriate to the client and used appropriately in a safe environment.

This is not a toy.

Use all appropriate safety precautions e.g. mats of an appropriate thickness, helmets

Check your suspension system regularly

This product is not a toy. Any activity involving motion, height, or rotation may cause serious injury. This product is intended for use under the supervision of trained professionals.

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