



MOVIVIT® MATS 20451, 20452, 20453



CONTENTS

1 x Movivit® Mat

USE

1. The mat can be used on either side.
2. Use on a flat, stable and appropriate surface. Use away from obstacles and hazards.

SAFETY

1. Always supervise the use of this product.

CARE AND MAINTENANCE

1. Store flat or roll up.
2. When stored, do not place heavy objects on the top.
3. Clean in lukewarm water with a mild cleaning agent. Do not use bleach or abrasive cleaning products.
4. Keep away from sources of heat and extreme temperatures.

SPECIFICATION

Closed cell foam

Fungicidal

Anti-bacterial

Stain-, sweat-, saliva- and water-resistant

Phthalate-free

Latex-free

Smooth surface

	Weight	Size
20451:	1.9kg	140 x 60 x 1.5cm
20452:	2.4kg	180 x 60 x 1.5cm
20453:	4.1kg	180 x 100 x 1.5cm

Further copies of these instructions can be downloaded at www.rompa.com

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at producthelp@rompa.com



COLCHÕES MOVIVIT 20451, 20452, 20453



CONTEÚDO

1 x Colchão Movivit

UTILIZAÇÃO

1. Pode utilizar ambas as faces do colchão.
2. Utilize numa superfície apropriada plana e estável. Utilize afastado de obstáculos e perigos.

SEGURANÇA

1. Guarde o colchão plano ou enrolado.
2. Quando estiver guardado, não coloque objetos pesados em cima.
3. Limpe com água morna e detergente macio. Não utilize lixívia ou produtos de limpeza abrasivos.
4. Mantenha afastado de fontes de calor e temperaturas extremas.

ESPECIFICAÇÕES

Células de espuma fechadas

Fungicida

Anti bacteriano

Resistente à água, manchas, doces e saliva

Sem ftalato

Sem látex

Superfície macia

	Peso	Tamanho
20451:	1.9kg	140 x 60 x 1.5cm
20452:	2.4kg	180 x 60 x 1.5cm
20453:	4.1kg	180 x 100 x 1.5cm